## Kaltlichtquellen

KL5110 / KL5120 / KL5125

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## Fiber-Optic Light Sources

KL5110 / KL5120 / KL5125

Operating Manual

## Inhalt

### KL5110

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## Index

### KL5110

1. Components	17
2. Assembly	17
3. Exchanging bulb	18
4. Maintenance	19
5. Technical data	19

### KL5120

1. Components	20
2. Assembly	20
3. Exchanging bulb	21
4. Maintenance	21
5. Technical data	22

### KL5125

1. Components	23
2. Assembly	23
3. Exchanging bulb	24
4. Maintenance	24
5. Technical data	25
Warranty conditions	26
Warranty extension	27

## KL5110

### 1. Components

1. Case
2. Bulb-holding device
3. Fiber-optic light guide  
(not included in offer)
4. Ventilator
5. Mains-supply socket
6. Transformer

### 2. Assembly

- Insert the power plug into the socket.
- Insert the transformer into the power source.
- Insert the fiber-optic light guide into the bulb-holding device.

### 3. Exchanging bulb

**Important advice!**

Disconnect the mains plug before every bulb exchange!

Fingerprints on the glass bulb impair the light efficiency and reduce the life of the bulb!

- Loosen the two screws and pull out the bulb-holding device.
- Take off the bulb socket.
- Loosen the bayonet joint.
- Remove the faulty bulb.  
Insert the new bulb between the iris diaphragm and the spring. Fix it with the bayonet joint. Pin up the bulb socket. Insert the bulb-holding device and fix it with the two screws.

### 4. Maintenance

- The instrument should only be used in a dry and well ventilated room. It should never come in contact with water or dampness.
- Do not drop the cold light source. High-precision electronic parts are extremely sensitive and may get broken if dropped. Also mechanical parts are made with highest precision and obviously will be damaged if dropped on the floor.
- In the event of a malfunction, your instrument should be repaired by A.KRÜSS Optronic.

### 5. Technical Data

Power supply primary voltage	100 – 240 V, 24 W
Power supply secondary voltage	7.5 V, 20 W
Halogen bulb	8 V, 20 W
Dimensions	170 x 100 x 55 mm
Weight	500 g

## KL5120

### 1. Components

1. Case
2. Bulb-holding device
3. Fiber-optic light guide (not included in offer)
4. Ventilator
5. Illumination socket
6. Mains switch on/off
7. Mains-supply socket
8. Power-supply unit

### 2. Assembly

- Remove instrument from package.
- Insert the power plug into the mains-supply socket and into the power source.
- Insert the fiber-optic light guide into the bulb-holding device.
- Adjust light intensity with the built-in iris diaphragm by turning the knob.

### 3. Exchanging bulb

#### *Important advice!*

Disconnect the mains plug before every bulb exchange!

Fingerprints on the glass bulb impair the light efficiency and reduce the life of the bulb!

- Loosen the three screws.
- Pull out the bulb-holding device.
- Loosen the bayonet joint.
- Remove the faulty bulb.  
Insert the new bulb between the iris diaphragm and the spring. Fix it with the bayonet joint.  
Pin up the bulb socket. Insert the bulb-holding device and fix it with the three screws.

### 4. Maintenance

- The instrument should only be used in a dry and well-ventilated room. It should never come to contact with water or dampness.
- Do not drop the cold light source. High-precision electronic parts are extremely sensitive and may get broken if dropped. Also mechanical parts are made with highest precision and obviously will be damaged if dropped on the floor.
- In the event of a malfunction, your instrument should be repaired by A.KRÜSS Optronic.

## 5. Technical Data

Power supply primary voltage	100 – 240 V, 20 W
Power supply secondary voltage	7.5 V, 20 W
Halogen bulb	8 V, 20 W
Dimensions	200 x 118 x 186 mm
Weight	2000 g
Specials	Intensity control without colour change Built-in adapter for spectroscope wavelength scale illumination

## KL5125

### 1. Components

1. Case
2. Bulb-holding device
3. Brightness control
4. Fiber-optic light guide (not included in offer)
5. Ventilator (at the bottom of the device)
6. Mains-supply socket
7. Fuse holder
8. Mains switch on/off

### 2. Assembly

- Remove instrument from package and set up the stands.

**Important notice!**  
Always set up the stands! The ventilator is located at the bottom of the device.

- Insert the power plug into the mains supply socket and into the power source.
- Insert the fiber-optic light guide into the bulb-holding device.
- Adjust the light intensity with the built-in iris diaphragm. Turn the knob to the left or right.

### 3. Exchanging the bulb

*Important advice!*

Disconnect the mail plug before every bulb exchange!

Fingerprints on the glass bulb impair the light efficiency and reduce the life of the bulb!

- Loosen the four screws.
- Remove the bulb-holding device.
- Remove the faulty bulb.  
Insert the new bulb between the cover and the spring.  
Insert the bulb-holding device and fix it with the four screws.

### 4. Maintenance

- The instrument should only be used in a dry and well-ventilated room. It should never come into contact with water or dampness. The instrument is equipped with a three-wire plug. This plug should only be used in sockets that are earthed.
- Do not drop the cold light source. High-precision electronic parts are extremely sensitive and may get broken if dropped. Also mechanical parts are made with highest precision and obviously will be damaged if dropped on the floor.
- In the event of a malfunction, your instrument should be repaired by A.KRÜSS Optronic.

### 5. Technical data

<b>Working voltage</b>	230 V, 150 W
<b>Halogen bulb</b>	15 V, 150 W
<b>Dimensions</b>	245 x 200 x 150 mm
<b>Weight</b>	4000 g

## Warranty conditions

A.KRÜSS Optronic guarantees the materials and workmanship of KRÜSS fiber-optic light sources for a period of 24 months from the date of shipping. A.KRÜSS Optronic will repair or replace defective devices within this period, given that they fall under terms of the guarantee. The device must be sent back to A.KRÜSS Optronic for warranty repairs or service. Shipping from the customer is at the expense of A.KRÜSS Optronic for warranty repairs. In all other cases, expenses must be paid by the customer.

A.KRÜSS Optronic guarantees that the hardware specified by A.KRÜSS Optronic for use with this device will function without error if used according to the instructions provided in our manual.

A.KRÜSS Optronic does not guarantee error-free and uninterrupted operation of the device or the accuracy of this instructions manual. A.KRÜSS Optronic is not liable for consequential damage.

### Warranty limitations:

This warranty does not cover errors and damage arising due to improper handling, through modification, misuse, operation above and beyond the specified environmental conditions or through unauthorised maintenance.

Further claims will not be accepted or recognised. A.KRÜSS Optronic expressly provides no guarantee for the workability or economic use in specific application cases.

A.KRÜSS Optronic reserves the right to make changes to this instruction manual and the technical data of the device in question at any point in time.

All devices may only be transported when properly packed in the complete original packaging. Request replacement packaging from your supplier if necessary.

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A.KRÜSS Optronic GmbH