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

**1. Stand ST 1512**

- 1 Polarization filter
- 2 Rotable glass plate
- 3 Conoscope lens
- 4 Lock lever
- 5 Rod
- 6 Spectroscope holder
- 7 Rotable disc for stone
- 8 Base

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1. Fix the spectroscope holder on the rod (5) with the screw.

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2. Fix the spectroscope with the screw.

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- Place a gem on  
the disc.  
Illuminate with the  
fiber light pipe.

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Note:

The spectroscope and the  
fiber light pipe have to form a  
90° angle.

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- Place the rotatable  
glass plate on  
the base.

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- Fix the conoscope (3)  
on the rod (5) with  
the screw.

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
- Fix the polarizing filter  
(1) on the rod (5) with  
the screw.

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- Insert the fiber  
light pipe into the side inlet  
to obtain an incident  
illumination.

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1. Hold the spectroscope with the slit side (1) in front of the medium to be examined.

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3. Focus the spectrum by sliding the slit tube (3).

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4. Adjust the slit width by turning the slit ring (2).

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5. The scale (7) is projected into the field of view by a 90° prism.

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6. Focus the scale by sliding the slit tube (6).

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**Applications:**

For teaching: Sun light reveals the solar spectrum known as “ Fraunhofers lines“  
 Light from a lamp shows the spectrum without lines


Absorption spectrum: Fluorescent light reveals the emission spectrum accordingly chromaticity diagram.

Emission spectrum: Sodium lamp shows spectral line at 589.3 nm

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
3 Fan

4 Mains supply socket

5 Spectroscope illumination socket

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**Assembly**

1. Remove the instrument from the package.
2. Insert the power plug into the mains supply socket (4) and the power socket.
3. Insert the flexible fiber light pipe into the lamp holding device (2).

4. Adjust the light intensity with the built-in iris diaphragm. Turn the knob to the right or left accordingly.

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
**Exchanging lamp****Disconnect the mains plug before every lamp exchange.**

Fingerprints on the glass lamp impair the light efficiency and reduce the life of the lamp.

1. Loosen the three screws. Pull out the bulb holding device.

2. Remove the lamp socket.

3. Loosen the bayonet joint

4. Remove the faulty lamp. Insert the new lamp between the iris diaphragm and the spring. Fasten it with the bayonet joint. Pin in the lamp socket. Insert the lamp-holding device and fasten it with the three screws.

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**Technical data**

Mains supply: 115-230V, 20W  
Halogen lamp: 8V, 20W  
Dimensions: 200x118x186 mm  
Weight: 2 kg

Intensity control without color change

8 V socket for spectroscopy wavelength-scale illumination

**3. Maintenance and care**

1. The instrument should only be used in a dry and well ventilated room. It should never come into contact with water or dampness. The instrument is equipped with a three wire plug. This plug should only be used in sockets that are earthed.
2. Do not drop the instruments. Optical and electronic parts are extremely sensitive and may get broken by dropping. Although mechanical parts are made with the highest precision they could be damaged if dropped.
3. Never touch the lens. Fingerprints are very difficult to remove. However, if they occur, use soft linen moistened with isopropyl alcohol to remove it.
4. In the event of a malfunction, your instrument should only be repaired by KRÜSS technicians.